

Quantum Healthcare
A ParagonCare Brand

CUREjet

TO MEET THE INCREASING DEMANDS FOR TARGETED PROCEDURES
BY MINIMIZING PAIN, ALLOW FASTER TREATMENT TIME,
AND PROVIDE PRECISE DELIVERY WITH MINIMUM DRUG LOSS

fast and precise, with less pain

BAZ BIOMEDIC'S CUREJET PROVIDES HIGH-SPEED CONSECUTIVE
& PRECISE NEEDLE-FREE INJECTION
ALLOWING REDUCED TREATMENT TIME AND TARGETS SPECIFIC APPLICATIONS
*ACNE SCAR, SCAR, STRETCH MARK, SCALP, MESOTHERAPY, ETC.



***Depth
control***

SELECTIVE DEPTH DELIVERY
EPIDERMIS OR DERMIS LAYERS



***Volume
control***

0.8 μ l ~ 1 μ l VOLUME
PER SHOT




***Less
pain***

200 μ m NOZZLE


PDLLA & HA INJECTION by CUREJET

BY CONTROLLING THE DEPTH OF DELIVERY, CUREJET ALLOWS THE DRUG SOLUTION TO PENETRATE THROUGH THE STRATUM CORNEUM OR THE DERMIS (0.05~2.5MM)


* 5mm Double Shots (FACE, H&E Stain x40)




* 5mm Double Shots (FACE, Alcian Blue Stain x40)




* 15mm 100 Shots (FACE, H&E Stain x40)



* 15mm 100 Shots (FACE, Alcian Blue Stain x40)



COMPATIBLE WITH VARIOUS DRUGS



CUREjet

POWER SUPPLY	220VAC, 60HZ, 200W
OUTPUT VOLTAGE	100~200V
REPETITION RATE	1, 3, 5, 10, 15, 20 HZ
DEMENSION	400[W] x 300[D] x 900[H]
WEIGHT	30kg

baz

바즈바이오메딕(주) BAZBIOMEDIC CO., LTD.

서울특별시 관악구 관악로 1, 서울대학교 연구공원 본관 (940동) 115호~119호
Room 115~119, Bldg 940, Seoul National Uni. 1, Gwanak-ro, Gwanak-gu, Seoul

www.bazbiomedic.com +82.2.883.1507